methylSig: A package for whole genome DNA methylation analysis

Yongseok Park, Maria E. Figueroa, Laura S. Rozek, and Maureen A. Sartor 2015-10-08

Contents

1	Intr	roduction	2
2	Inst	tallation	2
3	Basic usage		2
	3.1	Methylation score files	2
	3.2	Reading methylation score files	3
	3.3	Differential methylation analysis	4
	3.4	Site specific analysis	4
	3.5	Tiled data analysis	4
	3.6	Variance from one group	5
	3.7	Using local information	5
	3.8	Multi-thread computation	6
4	1 22	notation	6
4	4.1	CpG islands	6
		•	
	4.2	RefGene annotation	6
	4.3	Transcription factor (TF) enrichment test	8
	4.4	Data visualization	11
5	Dat	ta classes	12
	5.1	methylSigData object	12
		5.1.1 S4 data structure	12
		5.1.2 Subsetting	14
		5.1.3 Getting values	15
	5.2	methylSigDiff object	15
		5.2.1 S4 data structure	15
		5.2.2 Subsetting	15
		5.2.3 Getting values	16
		5.2.4 How to subtract DMCs or DMRs	16
	5.3	Summarizing data	16
	5 4		17

1 Introduction

DNA methylation plays critical roles in gene regulation and cellular specification without altering DNA sequences. It is one of the best understood and most intensively studied epigenetic marks in mammalian cells. Treatment of DNA with sodium bisulfite deaminates unmethylated cytosines to uracil while methylated cytosines are resistant to this conversion thus allowing for the discrimination between methylated and unmethylated CpG sites. Sodium bisulfite pre-treatment of DNA coupled with next-generation sequencing has allowed DNA methylation to be studied quantitatively and genome-wide at single cytosine site resolution.

methylSig is a method for testing for differential methylated cytosines (DMCs) or regions (DMRs) in wholegenome bisulfite sequencing (bis-seq) or reduced representation bisulfite sequencing (RRBS) experiments. methylSig uses a beta binomial model to test for significant differences between groups of samples. Several options exist for either site-specific or sliding window tests, combining strands, and for variance estimation. It allows annotating the resulting regions to multiple genome features, and visualizing the results for chosen genomic regions.

This document provides a step by step guide for the methylSig package.

2 Installation

methylSig is available on GitHub at http://www.github.com/sartorlab/methylSig, and the easiest way to install it is as follows:

```
library(devtools)
install_github('sartorlab/methylSig')
```

3 Basic usage

3.1 Methylation score files

methylSig expects input data to be formatted as follows:

```
##
           chrBase
                              base strand coverage freqC freqT
                      chr
    chr21.43008527 chr21 43008527
##
                                        F
                                                 32 100.00 0.00
                                        F
    chr21.43008531 chr21 43008531
                                                    96.88
##
    chr21.43008543 chr21 43008543
                                        F
                                                     90.62
                                                            9.38
    chr21.43008674 chr21 43008674
                                        R
                                                 27 100.00
                                                            0.00
##
##
    chr21.43008710 chr21 43008710
                                        R
                                                     94.03
                                                            5.97
    chr21.43008720 chr21 43008720
                                        R
                                                     92.54
```

Such CpG methylation score files can be obtained using bismark with the --bedGraph --cytosine_report flags in the bismark_methylation_extractor. The following awk command can be used on each cytosine report to obtain the correct format. Note, in \$4 + \$5 > 10, the 10 refers to the desired minimum coverage at each CpG site.

```
awk -v OFS="\t" '$4 + $5 > 10 {print $1"."$2, $1, $2, $3, $4 + $5, $4, $5}' in > out
```

The CpG methylation score file must contain at least seven columns. Among these, second to seventh column must be, in order, chromosome, base, strand, coverage, percentage of Cs and percentage of Ts. Column names are not important. Strand format is F/R or +/-, where F/+ represents forward and R/- represents reverse strand.

3.2 Reading methylation score files

methylSig package provides the methylSigReadData() function to read CpG methylation score files and convert these files into a methylSigData object for further analysis and annotation. The parameters and default options are:

```
methylSigReadData(fileList, sample.ids, assembly = NA, pipeline = NA,
  header = TRUE, context = NA, resolution = "base", treatment,
  destranded = TRUE, maxCount = 500, minCount = 10, filterSNPs = FALSE,
  num.cores = 1, quiet = FALSE)
}
```

Using data built into methylSig a typical read call might look like:

```
## Reading file (1/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/AML_1.txt
## (1/8) Count Invalid List: 0/2411=0
## Reading file (2/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/AML_2.txt
## (2/8) Count Invalid List: 11/4259=0.00258
## Reading file (3/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/AML_3.txt
## (3/8) Count Invalid List: 39/3861=0.0101
## Reading file (4/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/AML_4.txt
## (4/8) Count Invalid List: 2/3750=0.000533
## Reading file (5/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/NBM_1.txt
## (5/8) Count Invalid List: 0/6228=0
## Reading file (6/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/NBM_2.txt
## (6/8) Count Invalid List: 0/6430=0
## Reading file (7/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/NBM_3.txt
## (7/8) Count Invalid List: 1/5962=0.000168
## Reading file (8/8) -- /Users/raymond/Library/R/3.2/library/methylSig/extdata/NBM_4.txt
## (8/8) Count Invalid List: 23/3694=0.00623
## (1)(2)(3)(4)(5)(6)(7)(8)
```

It is possible for the user to filter out CpG sites based on the read coverage. CpG sites with very large read coverage may be due to PCR bias and hence including CpG sites with very high coverage may distort the statistics of data analysis. The methylSigReadData() function provides minCount and maxCount arguments for defining lower and upper limits for coverage. The default values are 10 and 500 respectively. It is also

possible to exclude C > T SNPs determined by the 1000 Genomes Project with the filterSNPs option. This is not done by default.

There are many arguments for the methylSigReadData() funciton. Among these fileList, sample.ids and treatment are required. Some options have default values, for example, destranded=TRUE, num.cores=1, and quiet=FALSE. Other arguments such as assembly, context and pipeline are optional and for information purposes only. The data type of treatment is a numeric vector. Each number represents a group. Multiple groups can be stored in one methylSigData object.

The arguemnt num.cores is used for multi-thread reading.

3.3 Differential methylation analysis

The main function of this package is the differential methylation analysis function methySigCalc(). It calculates differential methylation statistics between two groups of samples. It uses a beta-binomial approach to calculate differential methylation statistics, accounting for coverage and variation among samples within each group.

```
methylSigCalc(meth, groups = c(Treatment = 1, Control = 0),
   dispersion = "both", local.disp = FALSE, winsize.disp = 200,
   local.meth = FALSE, winsize.meth = 200, min.per.group = c(3, 3),
   weightFunc = methylSig_weightFunc, T.approx = TRUE, num.cores = 1)
}
```

3.4 Site specific analysis

The default is to do site specific analysis and to use both groups to estimate variances.

```
myDiffSigboth = methylSigCalc(meth, groups=c(1,0), min.per.group=3)
## Total number of bases: 3.26k
```

The differentially methylated cytosines (DMCs) can be defined based on qualues, pvalues and the methylation rate difference between two tested groups.

3.5 Tiled data analysis

methylSig package also provides methylSigTile() function to tile data within continuous non-overlapping windows. The default window size is 25bp. After tiling data, the methylSigCalc() function can be used to calculate differential methylation statistics.

```
### Tiled analysis
methTile = methylSigTile(meth,win.size = 25)
myDiffSigbothTile = methylSigCalc(methTile, groups=c(1,0), min.per.group=3)
```

```
## Total number of regions: 994
```

3.6 Variance from one group

Using the dispersion argument, it is possible to estimate variances from one group rather than from both groups. The following code calculates differential methylation statistics based on estimating variances from group 0 only.

3.7 Using local information

Total number of regions: 994

It is also possible to use information from nearby CpG sites to improve the variance and methylation level estimates. The default winsize.disp and winsize.meth are 200 bps. The winsize.disp argument only takes into effect when local.disp is set to TRUE'. Similarly \verb@winsize.meth@ argument only takes into effect when \verb@local.meth@ is set to TRUE'.

```
### Variance from both groups and using local information for variance
myDiffSigBothLoc = methylSigCalc(meth, groups=c(1,0),
         min.per.group=3, local.disp=TRUE, winsize.disp=200)
## Total number of bases: 3.26k
### Variance from sample treatment group "0" only and using local information for variance
myDiffSignormLoc = methylSigCalc(meth, groups=c(1,0), dispersion=0,
         min.per.group=3, local.disp=TRUE, winsize.disp=200)
## Total number of bases: 3.26k
### Variance from both groups and using local information for methylation level
myDiffSigBothMLoc = methylSigCalc(meth, groups=c(1,0),
         min.per.group=3, local.meth=TRUE, winsize.meth=200)
## Total number of bases: 3.26k
### Variance from both groups and using local information for methylation level and variance
myDiffSigBothMDLoc = methylSigCalc(meth, groups=c(1,0),
         min.per.group=3, local.disp=TRUE, winsize.disp=200,
         local.meth=TRUE, winsize.meth=200)
```

Total number of bases: 3.26k

3.8 Multi-thread computation

methylSig provides multicore programming to substantially reduce data analysis time. In the functions methylSigReadData and methylSigCalc, multi-core programming will be initiated using num.cores argument. Note that this option depends on R package parallel and hence is not available in the Windows platform. The following example illustrates the use of 2 cores.

4 Annotation

4.1 CpG islands

There are two functions, cpgAnnotation() and cpgAnnotationPlot(), in the methylSig package for CpG island annotation. The CpG island information file can be download the UCSC genome browser. The appropriate genome assembly should be used.

In Linux, the user may use the following command to download the annotation file for hg19. Please use appropriate directories for hg18, mm9 or mm10.

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/cpgIslandExt.txt.gz
gunzip *.gz
```

Here we use the CpG island annotation file provided in the methylSig package to annotate our example. Note that this is a reduced annotation file and is not appropriate for a full real data analysis.

```
cpgAnnotationPlot(cpgAnnDmc,main="DMCs")
```

4.2 RefGene annotation

Again, there are two functions, refGeneAnnotation() and refGeneAnnotationPlot(), in methylSig package for annotation using RefGene models. The refGene information file can be download from websites such

ALL




Figure 1: CpG annotation plots

DMCs




Figure 2: CpG annotation plots

ALL




Figure 3: refGene annotation plots

UCSC genome browser. The appropriate genome assembly (the same genome assembly of the provided data) should be used.

In a linux server, the user may use the following command to download the annotation file for hg19. Please use appropriate directories for hg18, mm9 or mm10.

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/refGene.txt.gz
gunzip *.gz
```

We use refGene annotation file provided in the methylSig package to annotate in our example. Note that this is a reduced annotation file and is not appropriate for the a full real data analysis.

4.3 Transcription factor (TF) enrichment test

The functions getTFBSInfo() and methylSig.tfbsEnrichTest() are provided for reading the TFBS information file and implementing transcription factor enrichment test.

UCSC genome browser provides TFBS conserved track for hg18 and hg19. The following linux server shell command can be used to download these files:

DMC




Figure 4: refGene annotation plots

```
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/tfbsConsSites.txt.gz
wget ftp://hgdownload.cse.ucsc.edu/goldenPath/hg19/database/tfbsConsFactors.txt.gz
gunzip *.gz
```

Here, tfbsConsSites.txt is tracking information and can be used directory in function getTFBSInfo(). The explanation of variable names is listed in file tfbsConsFactors.txt.

Another TFBS track is from ENCODE for hg18, hg19 and mm9. However, the methylSig package cannot use this type of track directly. We provide ENCODE TFBS track files that suitable for methylSig package at http://sartorlab.ccmb.med.umich.edu/software.

To identify which TFs have significant level of hypermethylation or hypomethylation across their binding sites, which could indicate whether the TF is having a weaker or stronger regulatory effect, respectively, we first tile all reads from regions to which a particular TF is predicted to bind. We then apply our beta-binomial model to the data for each TF to identify TFs with hyper- or hypo-methylated binding sites.

To achieve this, we provide function methylSigTileTFBS() to tile all data corresponding to the same TF.

```
methTileTFs = methylSigTileTFBS(meth, tfbsInfo)
myDiffTFs = methylSigCalc(methTileTFs, groups=c(1,0))
```

```
## Total number of TFs: 126
```




Figure 5: TFBS Enrichment

4.4 Data visualization

MethylSig offers a unique two-tiered visualization of the methylation data depending on the zoom level. When the chromosome range is large (>1 million bp), the visualization function does not show individual sample data.

```
methylSigPlot(meth, "chr21", c(43000000, 43500000), groups=c(1,0),
    cpgInfo=cpgInfo,refGeneInfo=refGeneInfo,
    myDiff=myDiffSigboth,tfbsInfo=tfbsInfo,tfbsDense=F,sigQ=0.05)
```

```
## Warning in methylSigPlot(meth, "chr21", c(4.3e+07, 43500000), groups =
## c(1, : Range of the drawing area is two wide, please reduce the range!
```




Figure 6: Data visualization over a large range

For narrow regions where at most $500~\mathrm{CpG}$ sites have data reads, users can visualize sample-specific coverage levels and % methylation at each site, together with group averages, significance levels and a number of genomic annotations.




Figure 7: Data visualization over a narrow range

```
methylSigPlot(meth, "chr21", c(43800000, 43900000), groups=c(1,0),
    cpgInfo=cpgInfo, refGeneInfo=refGeneInfo,
    myDiff=myDiffSigboth,tfbsInfo=tfbsInfo,tfbsDense=F,sigQ=0.05)
```

5 Data classes

5.1 methylSigData object

5.1.1 S4 data structure

The methylSig package uses S4 object. The contents of methylSigData@ can be shown using the show() function in R or just type the object itself.

```
## methylSigData object with 7,571 rows
##
        chr
                           end strand coverage1 numCs1 numTs1 coverage2 numCs2
               start
## 1 chr21 43000564 43000564
                                                     NA
                                                            NA
                                                                              NA
     chr21 43000820 43000820
                                              NA
                                                     NA
                                                            NA
                                                                       NA
                                                                              NA
      chr21 43008527 43008527
                                              32
                                                     32
                                                             0
                                                                      124
                                                                             122
## 4 chr21 43008531 43008531
                                              32
                                                     31
                                                                      125
                                                                             125
                                                             1
## 5 chr21 43008543 43008543
                                              32
                                                     29
                                                             3
                                                                      125
                                                                             117
## 6 chr21 43008665 43008665
                                              NA
                                                                       NA
                                                                              NA
                                                     NA
                                                            NA
     chr21 43008673 43008673
                                              27
                                                     27
                                                             0
                                                                       29
                                                                              25
## 8 chr21 43008709 43008709
                                              67
                                                             4
                                                                              26
                                                     63
                                                                       31
## 9 chr21 43008719 43008719
                                              67
                                                     62
                                                             5
                                                                       33
                                                                              33
## 10 chr21 43014041 43014041
                                    +
                                              NA
                                                     NA
                                                            NA
                                                                      207
                                                                             202
##
      numTs2 coverage3 numCs3 numTs3 coverage4 numCs4 numTs4 coverage5 numCs5
## 1
          NA
                    NA
                                   NA
                                                     13
                                                             0
                            NA
                                              13
                                                                       NA
## 2
          NA
                    NA
                            NA
                                   NA
                                              NA
                                                            NA
                                                     NA
                                                                       NΑ
                                                                              NA
## 3
           2
                    48
                            44
                                    4
                                              NA
                                                     NA
                                                            NA
                                                                       95
                                                                              95
## 4
           0
                    48
                            48
                                    0
                                              NA
                                                     NA
                                                            NA
                                                                       94
                                                                              94
## 5
           8
                    48
                            44
                                    4
                                              NA
                                                     NA
                                                            NA
                                                                       94
                                                                              81
## 6
                    NA
          NA
                            NA
                                   NA
                                              NA
                                                     NA
                                                            NA
                                                                      148
                                                                             148
## 7
           4
                    NA
                            NA
                                   NA
                                              NA
                                                     NA
                                                            NA
                                                                      146
                                                                             146
## 8
           5
                            NA
                    NA
                                   NA
                                              NA
                                                     NA
                                                            NA
                                                                      149
                                                                             149
## 9
           0
                    NA
                            NA
                                   NA
                                              NA
                                                     NA
                                                            NA
                                                                      149
                                                                             149
## 10
           5
                    190
                           184
                                    6
                                             228
                                                    228
                                                             0
                                                                      136
                                                                             132
##
      numTs5 coverage6 numCs6 numTs6 coverage7 numCs7 numTs7 coverage8 numCs8
## 1
          NA
                    NA
                            NA
                                   NA
                                              NA
                                                     NA
                                                            NA
                                                                       NA
## 2
                                                             0
          NA
                    NA
                            NA
                                   NA
                                              11
                                                     11
                                                                       NA
                                                                              NA
## 3
                                             108
          0
                    93
                            93
                                    0
                                                    104
                                                             4
                                                                       NA
                                                                              NA
## 4
          0
                    93
                            93
                                    0
                                             104
                                                     95
                                                             9
                                                                       NA
                                                                              NA
## 5
          13
                    93
                            93
                                             105
                                                     96
                                                             9
                                                                       NA
                                                                              NA
                                    0
## 6
           0
                    178
                           178
                                    0
                                             162
                                                    156
                                                             6
                                                                       73
                                                                              73
## 7
                                             162
                                                                       73
                                                                              73
           0
                    178
                           156
                                   22
                                                    161
                                                             1
## 8
           0
                   184
                           184
                                    0
                                             163
                                                    156
                                                             7
                                                                       74
                                                                              74
## 9
           0
                    NA
                           NA
                                   NA
                                             165
                                                    158
                                                             7
                                                                       74
                                                                              74
## 10
           4
                    155
                           135
                                   20
                                             203
                                                    198
                                                             5
                                                                       61
                                                                              61
##
      numTs8
## 1
          NA
## 2
## 3
          NA
## 4
          NA
## 5
          NA
## 6
           0
## 7
           0
## 8
## 9
## sample.ids: AML1 AML2 AML3 AML4 NBM1 NBM2 NBM3 NBM4
## treatment: 1 1 1 1 0 0 0 0
## destranded: TRUE
## resolution: base
```

options: maxCount=500 & minCount=10 & filterSNPs=FALSE & assembly=hg18 & context=CpG

Here, NA means there was no data at this base location on the related sample.

5.1.2 Subsetting

Data can be subset using matrix style operations. Row represents base location and each column is a sample. Below is an example to obtain data for samples 1 to 4:

```
meth1_4 = meth[,1:4]
```

This example returns the first 100 methylation reads in the data:

```
methSub1_100 = meth[1:100,]
```

Two arguments can be used together. This example returns the first 100 methylation reads for samples 1 and 2.

```
methSubData = meth[1:100,1:2]
methSubData
```

```
## methylSigData object with 60 rows
## ----
##
                           end strand coverage1 numCs1 numTs1 coverage2 numCs2
        chr
               start
## 1
     chr21 43008527 43008527
                                             32
                                                     32
                                                             0
                                                                      124
                                                                             122
## 2 chr21 43008531 43008531
                                             32
                                                             1
                                                                      125
                                                                             125
                                                     31
     chr21 43008543 43008543
                                             32
                                                     29
                                                             3
                                                                      125
                                                                             117
     chr21 43008673 43008673
                                             27
                                                             0
                                                                      29
## 4
                                                     27
                                                                              25
     chr21 43008709 43008709
                                             67
                                                     63
                                                             4
                                                                      31
                                                                              26
     chr21 43008719 43008719
                                             67
                                                     62
                                                                      33
                                                                              33
## 6
                                                             5
      chr21 43014041 43014041
                                             NA
                                                    NA
                                                                     207
                                                                             202
                                                            NA
     chr21 43014044 43014044
                                                            NA
                                                                     207
                                                                             195
                                             NA
                                                     NA
     chr21 43014076 43014076
                                             NA
                                                            NA
                                                                      202
                                                                             188
                                                     NA
## 10 chr21 43016439 43016439
                                             55
                                                     51
                                                                      93
                                                                              90
      numTs2
##
## 1
           2
## 2
           0
## 3
           8
## 4
           4
           5
## 5
## 6
           0
           5
## 7
## 8
          12
## 9
          14
## 10
           3
## sample.ids: AML1 AML2
## treatment: 1 1
## destranded: TRUE
## resolution: base
## options: maxCount=500 & minCount=10 & filterSNPs=FALSE & assembly=hg18 & context=CpG
```

5.1.3 Getting values

If the second argument is a string that matches one of the column names in the methylSigData object, it gives the values of that column. Valid column names are chr, start, end, strand, coverage1, ..., numCs1, ..., and numTs1.

```
coverage1 = meth[,"coverage1"]
startTop200 = meth[1:200,"start"]
```

5.2 methylSigDiff object

5.2.1 S4 data structure

The contents of methylSigDiff are:

```
myDiffSigboth
```

```
## methylSigDiff object with 3,260 rows
##
        chr
                          end strand
                                        pvalue
                                                  qvalue meth.diff
               start
## 1
     chr21 43008527 43008527
                                  + 0.30197133 0.6323745 -2.4167475
## 2
     chr21 43008531 43008531
                                  + 0.86910547 1.0000000 0.4158129
     chr21 43008543 43008543
## 3
                                  + 0.54717680 0.8571823 -2.8071020
## 4
     chr21 43014041 43014041
                                  + 0.35454483 0.6814532
                                                          2.3585424
## 5
     chr21 43014044 43014044
                                  + 0.91546697 1.0000000 -0.6133413
## 6
     chr21 43014076 43014076
                                  + 0.08082989 0.3398239 12.1681339
     chr21 43016439 43016439
## 7
                                  + 0.22017979 0.5537706 -2.1713344
## 8
     chr21 43016517 43016517
                                   - 0.38867776 0.7094566 -2.0930416
     chr21 43018213 43018213
                                  + 0.83266057 1.0000000 -1.4186090
## 10 chr21 43018274 43018274
                                   - 0.32935865 0.6514688 2.0009532
##
      logLikRatio
                        theta df
                                      mu1
## 1
      1.27484057 3.689198e+01 6 96.49842 98.91517
## 2
      0.02957797 1.172001e+01 6 98.45532 98.03950
      0.40674783 3.131946e+01 6 90.97520 93.78231
## 4
      0.98274395 2.829499e+01 7 97.82788 95.46934
## 5
      0.01210852 1.269081e+01 7 90.45700 91.07034
## 6
      4.39674527 1.904543e+01 6 93.23638 81.06825
## 7
      1.76894807 9.151486e+01 8 96.96983 99.14116
## 8
       0.86326517 1.000000e+06 6 93.10351 95.19656
## 9
       0.04765519 7.094155e+00 8 91.07918 92.49778
## 10 1.09884498 7.115526e+01 7 98.16913 96.16818
## -----
## sample.ids: AML1 AML2 AML3 AML4 NBM1 NBM2 NBM3 NBM4
## treatment: 1 1 1 1 0 0 0 0
## destranded: TRUE
## resolution: base
## options: dispersion=both & local.disp=FALSE & local.meth=FALSE& min.per.group=c(3,3)& Total: 3260
```

5.2.2 Subsetting

This object can also subset by row to obtain results from part of CpG sites or regions. However, the qualues will not be readjusted.

```
myDiff100 = myDiffSigboth[1:100,]
```

5.2.3 Getting values

Similar to the methSigData object, if the second argument is a string that is the same as one of the column names, it will return the results for that column. The valid variable names are chr, start, end, strand, pvalue, qvalue, meth.diff, logLikRatio, theta, df, mu1, and mu0. Here, for group methylation mean estimates mu1 and mu0, 1 and 0 come from the groups argument in the methylSigCalc() function. So if one has run the methylSigCalc() function with groups=c(4,0), then mu4 and mu0 will appear in the results.

```
qvalues = myDiffSigboth[,"qvalue"]
```

5.2.4 How to subtract DMCs or DMRs

This methylSigDiff object is very flexible to use by combining functions of subsetting and getting values. For example, the following code can obtain differentially methylated cytosines or regions defined as qvalue < 0.05 and difference of methylation rate > 25%.

Here abs() is the absolute value function in R.

If you want to use pvalues instead of qvalues, then you can use

5.3 Summarizing data

You can easily use other R functions to summarize or draw plots.

```
## (1/8) Count Invalid List: 0/2411=0

## (2/8) Count Invalid List: 0/4259=0

## (3/8) Count Invalid List: 0/3861=0

## (4/8) Count Invalid List: 0/3750=0

## (5/8) Count Invalid List: 0/6228=0

## (6/8) Count Invalid List: 0/6430=0

## (7/8) Count Invalid List: 0/5962=0

## (8/8) Count Invalid List: 0/3694=0

summary(methRaw[,"numCs1"]/methRaw[,"coverage1"])
```

```
## Min. 1st Qu. Median Mean 3rd Qu. Max. NA's
## 0.000 0.069 0.861 0.612 1.000 1.000 3606
```

```
summary(methRaw[,"coverage1"])
##
      Min. 1st Qu.
                    Median
                               Mean 3rd Qu.
                                               Max.
                                                       NA's
##
     10.00
             14.00
                     23.00
                              30.68
                                      39.00
                                             177.00
                                                        3606
hist(methRaw[,"numCs1"]/methRaw[,"coverage1"],
                 main="Histogram of methylation rate for sample 1",
                 xlab="methylation rate")
```

Histogram of methylation rate for sample




Figure 8: Methylation rate for sample 1

5.4 Generating heatmaps

Here we provide an example to generate a correlation heatmap.

Histogram of coverage for sample 1




Figure 9: Coverage for sample 1

Here is another example to generate a correlation heatmap based on differentially methylated cytosines.




Figure 10: Heatmap based on methylation rate at all CpG sites




Figure 11: Heatmap based on methylation rate at differentially methylated sites